 **Bramp**

# ICU TUBE FIXATION SYSTEM

---

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
# ICU TUBE FIXATION SYSTEM

## Indications for use

The ICU TUBE FIXATION SYSTEM (ICU TFS) has been developed for the fixation of single or dual lumen airway devices such as orotracheal tubes (OTT) and laryngeal masks (LM) or other devices (e.g. Combitube) after endotracheal or esophageal placement in order to reduce the risk of accidental extubation. ICU TFS is intended for this in individuals over the age of 8 years, for periods not exceeding 24 hours.

## ⚠ Warnings

A Warning states a condition, hazard or unsafe practice that could result in serious injury or death.

- This device should only be used by trained personnel in accordance with the User Guide.
- Apply just enough pressure against the tube to hold securely in place. Overtightening the elastic tie may restrict airflow.
- Do not use on children younger than 8 years of age.
- Designed for single patient use only.
- Do not reuse. Reuse will lead to increased risk of crosscontamination, degradation of performance and/or device malfunction. BRAMP is not responsible for any consequences of reuse.
- Before use, inspect all parts to ensure that they are free of damage and working properly.

## ⚠ Cautions

A Caution states a condition, hazard or unsafe practice that can result in minor personal injury or damage to the ICU TUBE FIXATION SYSTEM.

- The device should only be attached to the patient after insertion of an advanced airway device.
- Prior to the placement of the ICU TFS, the correct positioning of the OTT must be ensured.
- Do not use the device on a patient if it has been used during training.
- Do not disinfect or sterilize the product.
- The device is for short-term use only.
- Replace every 24h.
- Storage temperature: -34° C to +52° C.
- Operating temperature: -18° C to +43° C.

## Symbol Glossary

	According to the general safety and performance requirements of Regulation (EU) 2017/745 for medical devices
	Medical device
	Single patient use. Do not reuse
	Not sterile
	Not to be used in children younger than 8 years old
	Storage temperature (maximum and minimum)
	Warranty 2 years
	Keep dry
	Do not use if package is damaged
	Dispose of according to local protocol
	Not made with natural rubber latex

## Materials

Part	Material
Tube Anchoring Mask (TAM)	LDPE
Elastic Ties	TPS+SEBS
Interposition foam to TAM	MAYSER BLUEFOAM® DG-T + Double Coated Medical Polyester FilmTape
Fixation squares	70% Polyester - 30% Nylon
Other parts	Laminated Carded Mesh (45% PU – 39% PA – 16% PES)

## Notes:

The ICU TUBE FIXATION SYSTEM meets the following standards: EN 17989, ISO 14971.

The materials have not been tested for flammability; Non-sterile product.

Dispose of according to local protocol.

Information and Warranty: [www.bramp.com](http://www.bramp.com)

# Instructions for application and materials

## How to apply the Orotracheal Tube Fixation?

After placing the OTT (oro-tracheal tube), in the desired position in the tube anchoring mask (TAM) ③, place:


- a) The lower fixation strap ⑧ (below the occipital curvature);
- b) The tube anchoring mask (TAM) ③ is fixed with the fastening ends on both sides so that it is attached/fixed, but without putting too much pressure on the face;
- c) With the yellow elastic ties ④, fasten the OTT tube to the respective anchor pin of the TAM ③, pull firmly so that the OTT tube does not slide and turn so that it locks, leaving the elastic ties ④ under tension;

## How to fix the Gastric Probe?

d) The gastric probe can be fixed in various ways. It can be fastened to the pin that is free, or be fastened to the OTT, and so that there is no further doubt, a pink elastic tie ⑤ should be used for this purpose.

## How to Optimize the Fixation System?

- e) The following material are provided: a strap that should be placed perimetrally at the frontal-superior occipital level ⑦, with the function of keeping in position the inferior strap ⑧, being connected to it by the two occipital straps (fixed) and by two straps ⑥ each with 14 cm (loose straps) which are destined to the anterolateral fixation, that is, between the anterior temporal region above and laterally to the fixation point of the inferior strap to the TAM ③, inferiorly (in this way, the tendency for inferior mobilization of the OTT will be limited).
- f) A 25 cm strap ② is also supplied, which is intended to sagittally bind antero posteriorly the upper strap ⑦, ensuring that it does not slide downwards.


## Elements that compose the Kit

Fig	Part	Quantity x Designation
①		6x additional fixation squares for the sagittal and lateral accessory straps
②		1x 25 cm sagittal fixation accessory strap
③		1x tube anchoring mask (TAM)
④		4x yellow elastic ties
⑤		4x pink elastic ties
⑥		2x 14cm side facial fixation accessory straps
⑦		1x top fixation strap
⑧		1x lower fixation strap to TAM
⑨		1x interposition foam to TAM (this one goes glued to the TAM)

**Notes:** To reduce the risk of accidental extubation, it is strongly recommended that the colours are respected, thus, by convention, the yellow elastic ties hold the OTT and the pink elastic ties holds the N/OGT.

The use of any medical device is subject to the clinical conditions of the patient and as such subject to clinical judgement and good medical practice.

The figures are for illustrative purposes only, and may not be adjusted to the reality, or suffer slight changes.

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